

**HAZARDOUS STORAGE**  
**STOCK FLAMMABLE STORAGE CABINET**  
**88 SERIES**  
**88F944STOCK**



- High quality and reliable storage cabinet designed to maintain integrity for 30 minutes during fire
- Creates a safe and secure storage solution for flammable and hazardous materials when stored according to the material's own storage guidelines



Product Code	88F944STOCK
Hazard	Flammable
Standard Colour:	Golden Yellow (BS08E51)
H x W x D (mm):	915 x 459 x 459
Aperture Size (mm):	765 x 373
Weight (kg):	27
Steel Thickness (mm):	0.9
Standard No. of Shelves: Please note all shelves are perforated	1
Shelf Size (mm):	450 x 410 x 80
Additional Shelf Code:	88ES44FSTOCK
Optional Floorstand:	88FFS2 450 x 459 x 459
Number of Doors:	1— reinforced, rebated
Lock:	2 way locking chrome handle
Accessories Included:	Adjustable shelf Liquid tight sump tray
Shelf Capacity UDL (kg):	75
Certification:	DSEAR 2015 HS (G) 51 1990 Factory Inspectorate Certificate No 1 Parts 3&4
Hazard label:	Rec diamond, black border, black flame
Sump Tray Capacity (ltrs)	8.8
Sump Tray Size (mm):	360 x 410 x 80
Lead Time:	Immediate despatch



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

