

**HAZARDOUS STORAGE  
COSHH STORAGE CABINET  
88 SERIES  
88H844**



- High quality and reliable cabinet for the safe storage and segregation of substances hazardous to health based on employers carrying out mandatory COSHH risk assessment
- Can be used a precautionary measure to aid with the control and maintenance to harmful elements to help promote successful compliance with COSHH regulations 2002
- Substances should be stored within cabinets according to their own storage instructions in order to comply with COSHH regulations



Product Code	88H844
Hazard	COSHH
Standard Colour:	Silver Grey (7001)
H x W x D (mm):	1830 x 459 x 459
Aperture Size (mm):	1695 x 320
Weight (kg):	50
Steel Thickness (mm):	0.9
Standard No. of Shelves:	3
Shelf Size (mm):	450 x 385 x 40
Shelf Capacity (Itrs) :	6.5
Additional Shelf Code:	88ES44H
Number of Doors:	1— reinforced, rebated
Lock:	3 way locking chrome handle
Accessories Included:	Adjustable spillage tray Liquid tight sump tray
Shelf Capacity UDL (kg):	75
Certification:	DSEAR 2015: HS (G) 51 1990 Factory Inspectorate Certificate No 1 Parts 3&4
Hazard label:	White diamond, red border, black exclamation mark
Sump Tray Capacity (Itrs) :	12.5
Sump Tray Size (mm):	310 x 410 x 100
Lead Time:	15 days + delivery



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

**Wolflabs**

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

