


Cell and Tissue Culture Products



Cell and Tissue Culture Plates

Cell and Tissue Culture Plates are ideal for cell growth and cell yields on multiple, compare and other analyse.

- * Available with 8 different growth surface areas of 4, 6, 12, 24, 48, 96, 96U, 384 wells
- * Available in surface-treated or non-treated
- * Uniform well volume ensures equal growth surface area
- * Flat well bottom and round bottom plate
- * Well surface is smooth and free from striation to maximize usable growth area
- * Raised rims on wells with the uniform rings on the lid to reduce evaporation
- * Single position lid reduces the risks of cross-contamination and the handling mistakes
- * Wells are labeled with alphanumeric code for easy identification
- * Suitable for all common instruments and automation
- * Sterilized by gamma irradiation
- * DNase/RNase-free & Non-pyrogenic

Order Information


- * Individually packaged in peel-to-open paper/plastic blister pack
- * Every case has printed with lot No. for quality traceability



Cat.No	Well Qty	Surface Type	Max Well Volume(ml)[single well]	Lid	Sterile	Qty.per bag/case
TCP001004	4	General, Non-treated	1.50	Y	Y	1/100
TCP001006	6		17.0	Y	Y	1/100
TCP001012	12		6.80	Y	Y	1/100
TCP001024	24		3.50	Y	Y	1/100
TCP001048	48		1.55	Y	Y	1/100
TCP001096	96		0.39	Y	Y	1/100
TCP002096	96		0.32	Y	Y	1/100
TCP001384	384		0.14	Y	Y	1/100



Cat.No	Well Qty	Surface Type	Max Well Volume(ml)[single well]	Cell Growth Area(cm ²)[single well]	Lid	Sterile	Qty.per bag/case
TCP011004	4	Standard, surface treated	1.50	1.86-1.96	Y	Y	1/100
TCP011006	6		17.0	1.90-2.90	Y	Y	1/100
TCP011012	12		6.80	0.76-1.14	Y	Y	1/100
TCP011024	24		3.50	0.38-0.57	Y	Y	1/100
TCP011048	48		1.55	0.19-0.29	Y	Y	1/100
TCP011096	96		0.39	0.075-0.20	Y	Y	1/100
TCP012096	96U		0.32	0.32-0.38	Y	Y	1/100
TCP011384	384		0.14	0.10-0.11	Y	Y	1/100





Wolf Laboratories Limited

www.wolflabs.co.uk

Wolflabs

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

