


PIPETTE PRODUCTS

Pipets Products series, including serological pipets, aspirating pipets, transfer pipettes and rechargeable pipet aid.

All the pipets are manufactured with exclusively high-grade polystyrene(PS) or Polyethylene(PE) which are excellent for clear observation and reducing liquid attachment on the pipet surface to assure accurate delivery. The manufacturing environment is mastered by both ISO9001:2008, ISO13485 quality management and 100,000 grade clean-room cleanness system.

Furthermore, the universal design makes the pipets suitable for most brands of pipet-aid. Additionally, gamma irradiation sterilized of non-sterilized products are available in order to meet the different research needs and enhance the practicability. All the features guarantee them can be safely used in cell and tissue culture, bacteriological and clinical research.



Serological Pipets

- ★ Available with 7 capacities of 1.0, 2.0, 5.0, 10.0, 25.0 , 50.0ml and 100.0ml
- ★ 1.0, 2.0, 5.0ml are stretched, while 25.0, 50.0ml and 100.0ml are ultrasonically welded at tip and mouth piece
- ★ Both stretched mouth and ultrasonically welded mouth are available for 10ml pipet
- ★ Graduations are calibrated for accurate dispensing to within±2%
- ★ Different color rings are for easy identification
- ★ Bidirectional graduations on the pipets provide added applicability
- ★ Negative graduation allows additional working volume
- ★ All pipets are supplied with a filter plug
- ★ Strict consistency tested
- ★ Non-pyrogenic

Ordering Information

- ★ Individually packed in peel-to-open wrap with lot. No. for quality traceability
- ★ Carton has been specially designed for efficient bench top dispensing with a removable perforated front panel
- ★ Also available with bulk package in easy-to-open bag
- ★ Available with gamma irradiation sterilized or non-sterilized



Serological Pipets Individually Package (paper/plastic)

Cat.No	Volume[ml]	Length[mm]	Graduation[ml]	Color Code	Sterile	Package	Qty. Per bag/case
GSP010001	1	268.5±2.0	1/100	Yellow	Y	Paper/Plastic	500
GSP010002	2	272.0±2.0	1/50	Green	Y	Paper/Plastic	500
GSP010102	2	272.0±2.0	1/100	Green	Y	Paper/Plastic	500
GSP010005	5	341.0±2.0	1/10	Blue	Y	Paper/Plastic	200
GSP010010	10	346.3±2.0	1/10	Orange	Y	Paper/Plastic	200
GSP010110	10, Wide Mouth	346.3±2.0	1/10	Orange	Y	Paper/Plastic	200
GSP211010	10, Stretched Mouth	346.3±2.0	1/10	Orange	Y	Paper/Plastic	200
GSP010025	25	308.5±2.0	2/10	Red	Y	Paper/Plastic	150
GSP010125	25, Long	338.9±2.0	2/10	Red	Y	Paper/Plastic	150
GSP010050	50	346.2±2.0	5/10	Purple	Y	Paper/Plastic	100
GSP010100	100	346.8±2.0	1	Pink	Y	Paper/Plastic	50



Serological Pipets Individually Vacuum Package in bag (paper/plastic)

Cat.No	Volume(ml)	Length(mm)	Graduation(ml)	Color Code	Sterile	Package	Qty. Per bag/case
GSP110001	1	268.5±2.0	1/100	Yellow	Y	Paper/Plastic	100/600
GSP110002	2	272.0±2.0	1/50	Green	Y	Paper/Plastic	100/500
GSP110102	2	272.0±2.0	1/100	Green	Y	Paper/Plastic	100/500
GSP110005	5	341.0±2.0	1/10	Blue	Y	Paper/Plastic	50/200
GSP110010	10	346.3±2.0	1/10	Orange	Y	Paper/Plastic	50/200
GSP110110	10,Wide Mouth	346.3±2.0	1/10	Orange	Y	Paper/Plastic	50/200
GSP210010	10,Stretched Mouth	346.3±2.0	1/10	Orange	Y	Paper/Plastic	25/400
GSP110025	25	308.5±2.0	2/10	Red	Y	Paper/Plastic	50/150
GSP110125	25,Long	338.9±2.0	2/10	Red	Y	Paper/Plastic	50/150
GSP110050	50	346.2±2.0	5/10	Purple	Y	Paper/Plastic	30/90
GSP110100	100	346.8±2.0	1	Pink	Y	Paper/Plastic	10/50

Serological Pipets Bulk Vacuum Package

Cat.No	Volume(ml)	Length(mm)	Graduation(ml)	Color Code	Sterile	Qty. Per bag/case
GSP012001	1	268.5±2.0	1/100	Yellow	Y	25/1000
GSP012002	2	272.0±2.0	1/50	Green	Y	25/1000
GSP012102	2	272.0±2.0	1/100	Green	Y	25/1000
GSP012005	5	341.0±2.0	1/10	Blue	Y	25/500
GSP012010	10	346.3±2.0	1/10	Orange	Y	25/400
GSP012110	10, Wide Mouth	346.3±2.0	1/10	Orange	Y	25/400
GSP112010	10,Stretched Mouth	346.3±2.0	1/10	Orange	Y	25/400
GSP012025	25	308.5±2.0	2/10	Red	Y	10/150
GSP012125	25,Long	338.9±2.0	2/10	Red	Y	10/150
GSP012050	50	346.2±2.0	5/10	Purple	Y	10/100
GSP012100	100	346.8±2.0	1	Pink	Y	10/60
GSP011001	1	268.5±2.0	1/100	Yellow	N	25/1000
GSP011002	2	272.2±2.0	1/50	Green	N	25/1000
GSP011102	2	272.0±2.0	1/100	Green	N	25/1000
GSP011005	5	341.0±2.0	1/10	Blue	N	25/500
GSP011010	10	346.3±2.0	1/10	Orange	N	25/400
GSP011110	10, Wide Mouth	346.3±2.0	1/10	Orange	N	25/400
GSP111010	10,Stretched Mouth	346.3±2.0	1/10	Orange	N	25/400
GSP011025	25	308.5±2.0	2/10	Red	N	10/150
GSP011125	25, Long	338.9±2.0	2/10	Red	N	10/150
GSP011050	50	346.2±2.0	5/10	Purple	N	10/100
GSP011100	100	346.8±2.0	1	Pink	N	10/60

Serological Pipets Individually Package (plastic/plastic)

Cat.No	Volume(ml)	Length(mm)	Graduation(ml)	Color Code	Sterile	Package	Qty. Per case
GSP020001	1	268.5±2.0	1/100	Yellow	Y	Plastic/Plastic	500
GSP020002	2	272.0±2.0	1/50	Green	Y	Plastic/Plastic	500
GSP020102	2	272.0±2.0	1/100	Green	Y	Plastic/Plastic	500
GSP020005	5	341.0±2.0	1/10	Blue	Y	Plastic/Plastic	200
GSP020010	10	346.3±2.0	1/10	Orange	Y	Plastic/Plastic	200
GSP020110	10,Wide Mouth	346.3±2.0	1/10	Orange	Y	Plastic/Plastic	200
GSP021010	10,Stretched Mouth	346.3±2.0	1/10	Orange	Y	Plastic/Plastic	200
GSP020025	25	308.5±2.0	2/10	Red	Y	Plastic/Plastic	150
GSP020125	25,Long	338.9±2.0	2/10	Red	Y	Plastic/Plastic	150
GSP020050	50	346.2±2.0	5/10	Purple	Y	Plastic/Plastic	100
GSP020100	100	346.8±2.0	1	Pink	Y	Plastic/Plastic	50

Serological Pipets Individually Vacuum Package in bag (plastic/plastic)

Cat.No	Volume(ml)	Length(mm)	Graduation(ml)	Color Code	Sterile	Package	Qty. Per bag/case
GSP120001	1	268.5±2.0	1/100	Yellow	Y	Plastic/Plastic	100/600
GSP120002	2	272.0±2.0	1/50	Green	Y	Plastic/Plastic	100/500
GSP120102	2	272.0±2.0	1/100	Green	Y	Plastic/Plastic	100/500
GSP120005	5	341.0±2.0	1/10	Blue	Y	Plastic/Plastic	100/200
GSP120010	10	346.3±2.0	1/10	Orange	Y	Plastic/Plastic	50/200
GSP120110	10,Wide Mouth	346.3±2.0	1/10	Orange	Y	Plastic/Plastic	50/200
GSP121010	10,Stretched Mouth	346.3±2.0	1/10	Orange	Y	Plastic/Plastic	50/200
GSP120025	25	308.5±2.0	2/10	Red	Y	Plastic/Plastic	50/150
GSP120125	25,Long	338.9±2.0	2/10	Red	Y	Plastic/Plastic	50/150
GSP120050	50	346.2±2.0	5/10	Purple	Y	Plastic/Plastic	30/90
GSP120100	100	346.8±2.0	1	Pink	Y	Plastic/Plastic	10/50



- ★ Available with 7 capacities of 1.0, 2.0, 5.0, 10.0, 25.0, 50.0 and 100.0ml
- ★ 1.0, 2.0, 5.0, 10.0ml are stretched, while 10.0, 25.0, 50.0 and 100.0ml are ultrasonically welded at tip and mouthpiece
- ★ Colorless and transparent
- ★ Strict consistency tested
- ★ Non-pyrogenic

Ordering Information

- ★ Available with bulk package in easy to open plastic bags
- ★ Every case has printed lot No. for quality traceability
- ★ Available in sterilized or non-sterilized



Mini™ Serological Pipets

Mini pipets are specially designed for liquid handling operation in limited working space. These special designed pipets are significantly shorter than the normal pipet, only 9.25 inches. They come in three different volumes: 5.0ml, 10.0ml and 25.0ml.

Ordering Information

- ★ Available individually wrapped in an easy opening paper/plastic wrapped with lot number traceability
- ★ Inner boxes have been specially designed for efficient bench top dispensing with a removable perforated front panel
- ★ Also available with bulk package in easy to open bag

Cat.No	Volume(ml)	Graduation	Color Ring	Material	Sterile	Package	Qty. per bag/case
GSP010205	5.0	1/10	Blue	PS	Y	Paper/Plastic	1/200
GSP010210	10.0	2/10	Orange	PS	Y	Paper/Plastic	1/150
GSP010225	25.0	5/10	Red	PS	Y	Paper/Plastic	1/100

Aspirating Pipets

Aspirating pipets are ideal for liquid transfer of mix.

Cat.No	Capacity(ml)	Length(mm)	Material	Sterile	Package	Qty. per bag/case
GSP000001	1.0	268.5±2.0	PS	Y	Plastic Bag	25/1000
GSP000002	2.0	272.0±2.0	PS	Y	Plastic Bag	25/1000
GSP000005	5.0	341.0±2.0	PS	Y	Plastic Bag	25/500
GSP000010	10.0	346.3±2.0	PS	Y	Plastic Bag	25/400
GSP000025	25.0	308.5±2.0	PS	Y	Plastic Bag	10/150
GSP000125	25.0,Long	338.9±2.0	PS	Y	Plastic Bag	10/150
GSP000050	50.0	346.2±2.0	PS	Y	Plastic Bag	10/100
GSP000100	100	346.8±2.0	PS	Y	Plastic Bag	10/60
GSP001001	1.0	268.5±2.0	PS	N	Plastic Bag	25/1000
GSP001002	2.0	272.0±2.0	PS	N	Plastic Bag	25/1000
GSP001005	5.0	341.0±2.0	PS	N	Plastic Bag	25/500
GSP001010	10.0	346.3±2.0	PS	N	Plastic Bag	25/400
GSP001025	25.0	308.5±2.0	PS	N	Plastic Bag	10/150
GSP001125	25.0,Long	338.9±2.0	PS	N	Plastic Bag	10/150
GSP001050	50.0	346.2±2.0	PS	N	Plastic Bag	1/200
GSP001100	100	346.8±2.0	PS	N	Plastic Bag	1/200

10ml Aspirating pipets (stretched)

Cat.No	Capacity(ml)	Length(mm)	Material	Sterile	Package	Qty. per bag/case
GSP002010	10.0	346.3±2.0	PS	Y	Plastic	25/400
GSP003010	10.0	346.3±2.0	PS	N	Plastic	25/400
GSP101010	10.0	346.3±2.0	PS	Y	Paper/Plastic	50/200
GSP201010	10.0	346.3±2.0	PS	Y	Paper/Plastic	50/200



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.