


2xOFI 202-40M (oil-less)

Voltage	V	230	230	3x400	3x400														
Frequency	Hz	50	60	50	60														
Motor	HP	4.00	4.00	4.00	4.00														
	kW	2.94	2.94	2.94	2.94														
Displacement	l/min	580	658	580	658														
	CFM	20.48	23.24	20.48	23.24														
FAD @ 8 bar/120 psi	l/min	260	292	260	292														
	CFM	9.20	10.30	9.20	10.30														
Max pressure ¹⁾	bar	8	8	8	8														
	psi	120	120	120	120														
Max. current	A	16	18	10	10														
Tank size	litres	40	40	40	40														
	gallon	10.6	10.6	10.6	10.6														
Weight	kg	162	162	162	162														
	lbs	357	357	357	357														
Noise level	db(A)	63	66	63	66														
Thermal protection		Yes	Yes	Yes	Yes														
Duty cycle		100%	100%	100%	100%														



Dimensions (l x w x h)
661 x 774 x 856 mm / 26.0 x 30.5 x 33.7 inch

¹⁾ Higher pressure available upon request

Technical modifications reserved

JUN-AIR 

IDEX
IDEX CORPORATION





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

