


2xOFI202-90B (oil-less)

Voltage	V	200	200	230	240	3x400	3x440						
Frequency	Hz	50	60	50	60	50	60						
Motor	HP	4.00	4.00	4.00	4.00	4.00	4.00						
	kW	2.94	2.94	2.94	2.94	2.94	2.94						
Displacement	l/min	580	656	580	656	580	656						
	CFM	20.48	23.17	20.48	23.17	20.48	23.17						
FAD @ 8 bar ³⁾	l/min	260	292	260	292	260	292						
	CFM	9.18	10.31	9.18	10.31	9.18	10.31						
Max. pressure ¹⁾	bar	8	8	8	8	8	8						
	psi	120	120	120	120	120	120						
Max. current	A	-	-	16.0	18.0	10.0	10.0						
Tank size	litres	90	90	90	90	90	90						
	gallon	23,7	23,7	23,7	23,7	23,7	23,7						
Weight	kg	115	115	115	115	115	115						
	lbs	253	253	253	253	253	253						
Noise level	dB(A)/1m	79	78	79	78	79	78						
Thermal protection		Yes	Yes	Yes	Yes	Yes	Yes						
Duty cycle		100%	100%	100%	100%	100%	100%						



Dimensions (l x w x h)
1000 x 530 x 850 mm / 39.4 x 19.7 x 33.5 inch

- ¹⁾ Higher pressure available upon request
³⁾ Displacement is reduced by approx. 18-20% on units with dryer (D). Min. pressure required to operate dryer: 6 bar.

Technical modifications reserved



JUN-AIR 



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

