

## TECHNICAL FEATURE

4DL 31NBPR 00 705



### The :

- Seat and backrest slick surfaces
- Soft
- Easy to disinfect

### ADJUSTABLE CHAIRS / WIDE VINYL CHAIRS

Polyamide, 5-branch, reinforced base with heavy-duty nylon casters

Black gas-lift column with height adjustment from 470 to 600 mm

and shock absorber

Wide seat in grey M1 vinyl with high density soft foam

High backrest in grey M1 vinyl with high density soft foam,

adjustable in height and depth



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

