

## TECHNICAL FEATURE

4DL 40GBLP 00 705



### The :

- Seat and backrest slick surfaces
- Soft
- Easy to disinfect
- Multiple adjustments
- Fully adjustable

### ASYNCHRONOUS CHAIRS / WIDE VINYL CHAIRS

Chrome, 5-branch, scuff-resistant, reinforced base with articulated glides

Chrome gas-lift column with height adjustment from 410 to 540 mm

and shock absorber

Wide seat in grey M1 vinyl with high density soft foam

Independently adjustable seat and backrest tilt with one fingertip lever

thanks to an asynchronous mechanism

High backrest in grey vinyl with high density soft foam, adjustable in height



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

