

TECHNICAL FEATURE

4DL 40GBLR 00 705



The :

- Seat and backrest slick surfaces
- Soft
- Easy to disinfect
- Multiple adjustments
- Fully adjustable

ASYNCHRONOUS CHAIRS / WIDE VINYL CHAIRS

Chrome, 5-branch, scuff-resistant, reinforced base with heavy-duty nylon casters

Chrome gas-lift column with height adjustment from 450 to 580 mm

and shock absorber

Wide seat in grey M1 vinyl with high density soft foam

Independently adjustable seat and backrest tilt with one fingertip lever

thanks to an asynchronous mechanism

High backrest in grey vinyl with high density soft foam, adjustable in height



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

