

## TECHNICAL FEATURE

4GM 34GHLP 02 905



### The :

- Compact
- Polyvalent
- Easily adjustable

### STOOLS / VINYL COMFORT

Chrome, 5-branch, scuff-resistant, reinforced base with articulated glides  
Chrome gas-lift column with height adjustment from 590 to 830 mm and shock absorber  
Height adjustable polyamide foot-ring  
Round seat in black M1 vinyl with high density soft foam  
Articulated backrest in black M1 vinyl with high density soft foam, adjustable in height and depth



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

