


TECHNICAL FEATURE

4GM 37GPLR 00 905



The  :

- Compact
- Polyvalent
- Easily adjustable
- Hygienic

STOOLS / VINYL COMFORT

Chrome, 5-branch, scuff-resistant, reinforced base with heavy-duty nylon casters
Chrome gas-lift column with height adjustment from 600 to 700 mm by foot control and shock absorber
Round seat in black M1 vinyl with high density soft foam
Articulated backrest in black M1 vinyl with high density soft foam, adjustable in height and depth



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

