

## TECHNICAL FEATURE

4NN 10NHNR 02 905



### STOOLS / VINYL CLASSIC

Polyamide, 5-branch, reinforced base  
with heavy-duty nylon casters  
Black gas-lift column with height adjustment from 620  
to 860 mm  
and shock absorber  
Height adjustable polyamide foot-ring  
Round seat in black vinyl with high density soft foam

#### The + :

- Compact
- Polyvalent
- Easy to disinfect
- Easily adjustable



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

