

TECHNICAL FEATURE

7NG 35NHPP 02 905



The :

- Comfortable and airy surface
- Without staples nor seams
- Abrasion-resistant
- Easily adjustable

ADJUSTABLE CHAIRS / POLYURETHANE

Polyamide, 5-branch, reinforced base with articulated glides
Black gas-lift column with height adjustment from 580 to 820 mm and shock absorber
Height adjustable polyamide foot-ring
Ergonomic seat in black KANGO soft polyurethane
Articulated backrest in black soft polyurethane, adjustable in depth and height



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

