

TECHNICAL FEATURE

7NG 40GHLP 01 905



The :

- Comfortable and airy surface
- Without staples nor seams
- Abrasion-resistant
- Easily adjustable
- Fully adjustable

ASYNCHRONOUS CHAIRS / POLYURETHANE

Chrome, 5-branch, scuff-resistant, reinforced base with articulated glides
Chrome gas-lift column with height adjustment from 560 to 800 mm and shock absorber
Height adjustable aluminium foot-ring
Ergonomic seat in black KANGO soft polyurethane
Independently adjustable seat and backrest tilt with one fingertip lever
thanks to an asynchronous mechanism
Backrest in black soft polyurethane, adjustable in height



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

