



Low Temperature Chest Freezer CLT140

CLT140

- ✓ -10°C to -45°C operating temperature range
- ✓ External digital controller
- ✓ Lock and two 2 keys
- ✓ Spark free
- ✓ Audible and visual out of temperature alarms

SPECIFICATION

Dimensions	720mm W x 885mm H x 600mm D
Weight	59kg
Product name	Low Temperature Chest Freezer CLT140
GDPA code	CLT140
SKU code	444440160
External dimensions	720mm W x 885mm H x 600mm D
Internal dimensions	560mm W x 650mm H x 440mm D
Packed dimensions	960mm W x 960mm H x 700mm W
Product weight	59kg
Packed weight	59kg
Gross capacity	140L
Net capacity	133L
Colour	White
Door type	Solid
Door style	Hinged
Replacement door seal	N

Open door sensor	N
Lock	Y
Adjustable feet	N
Temperature range	-10°C to -45°C
Temperature display on panel	Y
Temperature alarm	Y
Interior light	N
Auto or manual defrost	Manual
Fan	N
Controller type	Digital, external
Climate class	N
Refrigerant type	R404
CFC free	Y
UK 3 pin plug	Y
Voltage supply	230/240
Spark free	Y
Auxiliary output	Y
Energy consumption	2.2 kW/24h



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

