



Low Temperature Chest Freezer CLT383

CLT383

- ✓ -10°C to -45°C operating temperature range
- ✓ External digital controller
- ✓ Lock and 2 keys
- ✓ Spark free
- ✓ Audible and visual out of temperature alarms

SPECIFICATION

Dimensions	1560mm W x 885mm H x 600mm D
Weight	80kg
Product name	Low Temperature Chest Freezer CLT383
GDPA code	CLT383
SKU code	444440162
External dimensions	1560mm W x 885mm H x 600mm D
Internal dimensions	1400mm W x 650mm H x 440mm D
Packed dimensions	1600mm W x 960mm H x 700mm D
Product weight	80kg
Packed weight	95kg
Gross capacity	383L
Net capacity	368L
Colour	White
Door type	Solid
Door style	Hinged
Replacement door seal	N

Open door sensor	N
Lock	Y
Adjustable feet	N
Temperature range	-10°C to -45°C
Temperature display on panel	Y
Temperature alarm	Y
Interior light	N
Auto or manual defrost	Manual
Fan	N
Controller type	Digital, external
Climate class	N
Refrigerant type	R404
CFC free	Y
UK 3 pin plug	Y
Voltage supply	230/240
Spark free	Y
Auxiliary output	Y
Energy consumption	5.2 kW/24h



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

