



LSC119UK ATEX Under-counter Laboratory Fridge-Freezer 119L

LM LSC119UK ATEX

- ✓ 2°C to 10°C and -18°C to -25°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms

SPECIFICATION

Dimensions	540mm W x 860mm H x 600mm D
Weight	34kg
Product name	Under-counter Laboratory Fridge-Freezer 119L LSC119UK ATEX
GDPA code	LM LSC119UK ATEX
SKU code	UK 444442777, EU 444441955
External dimensions	540mm W x 860mm H x 600mm D
Internal dimensions	Fridge 450mm W x 525mm H x 450mm D (step 225mm H x 155D), Freezer 400mm W x 120mm H x 295mm D
Packed dimensions	602mm W x 877mm H x 639mm D
Product weight	34kg
Packed weight	36kg
Gross capacity	119L (fridge 106L, freezer 13L)
Net capacity	117L (fridge 105L, freezer 12L)
External finish	White
Colour	White/green
Internal finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door style	Hinged, reversible
Door open depth	1090mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	2 front, 2 rear rollers
Temperature range	Fridge 2°C to 10°C, freezer -18°C to -25°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	2
Shelf material	Wire
Number of shelf positions	4

Shelf size	1x 308mm W x 4.5mm H x 448mm D, 1x 200mm W x 4.5mm H x 448mm D
Moulded drawer	1x solid
Adjustable door bins	2x full, 2x half
Interior light	N
Auto or manual defrost	Fridge auto, freezer manual
Fan	N
Climate class	N
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
UK 3 pin plug	Y
Voltage supply	230
Spark free	Y
Auxiliary output	Y
Insulation thickness	45mm
Insulation type	Cyclopentane
Air tempertaure alarm	Y
Stacking kit compatible	N





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

