



## UNDER-COUNTER LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 119L WHITE



- ✓ **Operating Temperature Range of +2° C to +10° C for Fridge and -18° C to -25° C for Freezer \***  
Conforms to MHRA guidelines
- ✓ **NEW Advanced Controller**  
Records temperatures accurate to 0.1° C
- ✓ **NEW Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**  
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**  
For ease of movement
- ✓ **Wire Shelves**  
Enables optimum air flow inside unit



[www.facebook.com/Lec-Medical](https://www.facebook.com/Lec-Medical)  
[www.twitter.com/Lec\\_Medical](https://www.twitter.com/Lec_Medical)

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440972	LSC119UK	Under-counter Laboratory Fridge-freezer Combi Solid Door 119L White

\* These are the set temperatures. Please contact us if you require guaranteed chamber temperatures.



## FULL SPECIFICATION – LSC119UK UNDER-COUNTER LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 119L WHITE



Description	Under-counter
SKU Code	444440972
EAN Code	5052263009721
GDPA Code	LSC119UK
Gross Capacity (L)	119 (Fridge: 106, Freezer: 13)
Net Capacity (L)	117
Unit External Width (mm)	540
Unit External Depth (mm)	600
Unit External Height (mm)	860
Product Weight (kg)	34
Packed Width (mm)	602
Packed Depth (mm)	639
Packed Height (mm)	877
Packed Weight (kg)	36
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	Fridge: 450 x 450 x 525, Freezer: 400 x 295 x 120 (step: 155d x 225h)
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	Fridge: +2° C to +10° C, Freezer: -18° C to -25° C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1090
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet	2 Adjustable Feet, 2 Rear Rollers
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	2
Number of Shelf Positions	4
Shelf Size	1: 308 x 448 x 4.5, 1: 200 x 448 x 4.5
Shelf Material	Wire
Moulded Drawers	1 Solid
Adjustable Door Bins	2 full + 2 half
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Fridge: Automatic, Freezer: Manual
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	45
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

