



## LSC324UK ATEX Freestanding Laboratory Fridge-Freezer 324L

### LSC324UK ATEX

- ✓ 2°C to 10°C and -18°C to -25°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms

## SPECIFICATION

Dimensions	600mm W x 1870mm H x 600mm D
Weight	70kg
Product name	Freestanding Laboratory Fridge-Freezer 324L LSC324UK ATEX
GDPA code	LSC324UK ATEX
SKU code	UK 444442890, EU 444441957
External dimensions	600mm W x 1870mm H x 600mm D
Internal dimensions	Fridge 510mm W x 805mm H x 430mm D (step 260mm H x 250mm D), Freezer 463mm W x 780mm H x 445mm D
Packed dimensions	657mm W x 1885mm H x 734mm D
Product weight	70kg
Packed weight	78kg
Gross capacity	324L (fridge 190L, freezer 134L)
Net capacity	302L
External finish	White
Colour	White/green
Internal finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door style	Hinged, reversible
Door open depth	1190mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	2 front, 2 rear rollers
Temperature range	Fridge 2°C to 10°C, freezer -18°C to -25°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	3
Shelf material	Wire
Number of shelf positions	6

Shelf size	320mm W x 6mm H x 505mm D
Moulded drawer	2x half solid
Adjustable door bins	2x full, 4x half
Interior light	N
Auto or manual defrost	Fridge auto, freezer manual
Fan	N
Climate class	N
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
UK 3 pin plug	Y
Voltage supply	230
Spark free	Y
Auxiliary output	Y
Insulation thickness	60mm
Insulation type	Cyclopentane
Air temperature alarm	Y
Stacking kit compatible	N





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

