



Under-Counter Laboratory Control Plus Fridge-Freezer Combi 119L

GDPA code: LSFC119UK Lab SF SD UC F/F

SKU code: 444410684

- ✓ Operating Temperature Range of +2°C to +10°C for Fridge and -18°C to -25°C for Freezer
- ✓ Spark Free
- ✓ Advanced Controller
- ✓ Unique Anti Microbial Handle
- ✓ Fitted Lock as Standard




The Experts in
Medical Refrigeration

Specification: Under-Counter Laboratory Control Plus Fridge-Freezer Combi 119L

| | |
|--|---|
| Gross Capacity (litres) | 119 (Fridge: 106, Freezer: 13) |
| Net Capacity (litres) | 117L (fridge 105L, freezer 12L) |
| Unit External Width (mm) | 540 |
| Unit External Depth (mm) | 600 |
| Unit External Height (mm) | 860 |
| Product Weight (kg) | 34 |
| Packed Width (mm) | 602 |
| Packed Depth (mm) | 639 |
| Packed Height (mm) | 877 |
| Packed Weight (kg) | 36 |
| Country of Origin | UK |
| Internal Dimensions (mm) (WxDxH) | Fridge: 450 x 450 x 525, Freezer: 400 x 295 x 120 (step: 155d x 225h) |
| Temperature Range | Fridge: +2°C to +10°C, Freezer: -18°C to -25°C |
| Colour | White/Green |
| External Finish | White |
| Internal Finish | White |
| Type of plastic | SRL800 (HIPS/ABS Mix) |
| Door Finish | Painted Steel / Solid |
| Door Type | Solid |
| Door Style | Hinged, Reversible |
| Door Open Depth (mm) | 1090 |
| Replacement Door Seal | No |
| Open Door Sensor | Yes |
| Lock | Yes |
| Adjustable Feet (number & position) | 2 Adjustable, 2 Rear Rollers |
| Porthole | No |
| Porthole Size | N/A |
| Porthole Location | N/A |
| Temperature Display on Panel | Yes |
| Temperature Alarm | Yes |
| Remote Alarm Contacts | Yes |
| Audible Alarms | Yes |
| Visual Alarms | Yes |
| Number of Shelves | 2 |
| Number of Shelf Positions | 4 |
| Shelf Size (mm) (WxDxH) | 308 x 448 x 4.5, 1: 200 x 448 x 4.5 |
| Shelf Material | Wire |
| Number of Baskets | 0 |
| Basket Material | 1 Solid |
| Basket Size (mm) (WxDxH) | 2 full + 2 half |
| Interior Light | N/A |
| Light Type | No |
| Light Activation | - |
| Fan | - |
| Spark Free | Yes |
| Auto or Manual Defrost | Yes |
| Controller Type | Fridge: Automatic, Freezer: Manual |
| Number of Probes | Digital, External |
| Probe Type | 1 |
| Climate Class | Air |
| Refrigerant Type | N |
| CFC Free | R600a |
| HFC Free | Yes |

Specification: Under-Counter Laboratory Control Plus Fridge-Freezer Combi 119L



| | |
|---|--------------|
| HCFC Free | Yes |
| Ammonia Free | Yes |
| UK 3 Pin Plug | Yes |
| Other Plug Types Available | Yes |
| Voltage Supply | - |
| Insulation Thickness (mm) | Yes |
| Insulation Type | 45 |
| Blowing Agent | Cyclopentane |
| Anti Microbial Handle | Yes |
| Optional Drawer System Available | No |
| Wall Brackets | No |
| Digital Lock | No |





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

