



UNDER-COUNTER LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 138L WHITE



- ✓ **Operating Temperature Range of +2° C to +10° C for Fridge and -12C to -20° C for Freezer**
Conforms to MHRA guidelines
- ✓ **Advanced Controller**
Records temperatures accurate to 0.1° C
- ✓ **Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**
For ease of movement
- ✓ **Wire Shelves**
Enables optimum air flow inside unit

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444411207	LSC138UK	Under-counter Laboratory Fridge-freezer Combi Solid Door 138L White



FULL SPECIFICATION – LSC138UK UNDER-COUNTER LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 138L WHITE

Description	Under-counter
SKU Code	444411207
EAN Code	5052263112070
GDPA Code	LSC138UK
Gross Capacity (L)	138 (Fridge: 120, Freezer: 18)
Net Capacity (L)	138 (Fridge: 120, Freezer: 18)
Unit External Width (mm)	595
Unit External Depth (mm)	660
Unit External Height (mm)	850
Product Weight (kg)	34
Packed Width (mm)	640
Packed Depth (mm)	660
Packed Height (mm)	910
Packed Weight (kg)	38
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	Fridge: 500 x 450 x 525, Freezer: 440 x 295 x 155
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	Fridge: +2° C to +10° C, Freezer: -12° C to -20° C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1130
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes (2 keys)
Adjustable Feet	4 Adjustable Feet
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	3
Number of Shelf Positions	5
Shelf Size	495 x 325 x 6
Shelf Material	Wire
No. of Baskets	1
Basket Size (mm) (W x D x H)	485 x 180 x 171
Adjustable Door Bins	3 full
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Fridge: Automatic, Freezer: Manual
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
UK Plug	Yes
Port Hole	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	45
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

