



PE109C Countertop Pharmacy Refrigerator Solid Door 45L

PE109C

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 2 year warranty as standard

SPECIFICATION

Dimensions	470mm W x 530mm H x 465mm D
Weight	20kg
Product name	Pharmacy Fridge Solid Door 45L PE109C
GDPA code	PE109C
SKU code	UK 444441788, EU 444442327,
External dimensions	470mm W x 530mm H x 465mm D
Internal dimensions	410mm W x 400mm H x 270mm D
Packed dimensions	554mm W x 565mm H x 586mm D
Product weight	20kg
Packed weight	23kg
Gross capacity	45L
Net capacity	37L
External finish	Painted steel
Colour	White/green
Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door finish	Painted steel
Door style	Hinged reversible
Door open depth	910mm
Replacement door seal	N
Open door sensor	N
Lock	Y
Adjustable feet	2 front
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	1
Shelf material	Wire
Number of shelf positions	1

Shelf size	402mm W x 6mm H x 200mm D
Interior light	N
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	N
Refrigerant type	R134a
CFC free	Y
HFC free	N
HCFC free	Y
Ammonia free	Y
Optional plug types	Y
Voltage supply	230
Spark free	N
Auxiliary output	N
Insulation thickness	30mm
Insulation type	Cyclopentane
Noise level	<45dB
Energy consumption	1.68 kw/24h
Door open alarm	Y
Air temperature alarm	Y
Antimicrobial handle	Y
Optional drawer system	N
Stacking kit compatible	N





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

