



## Countertop Pharmacy Refrigerator Solid Door 41L PE102C

PE102C

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 2 year warranty as standard

### SPECIFICATION

Dimensions	379mm W x 702mm H x 500mm D
Weight	23.5kg
Product name	Pharmacy Fridge Solid Door 41L PE102C
GDPA code	PE102C
SKU code	UK 444441789, EU 444442329, AU 444442330
External dimensions	379mm W x 702mm H x 500mm D
Internal dimensions	245mm W x 400mm H x 310mm D
Packed dimensions	450mm W x 740mm H x 560mm D
Product weight	23.5kg
Packed weight	25.5kg
Gross capacity	41L
Net capacity	35L
External finish	Painted steel
Colour	White/green

Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door finish	Painted steel
Door style	Hinged right hand only, non reversible
Door open depth	810mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	N
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	2
Shelf material	Wire
Number of shelf positions	5
Shelf size	270mm W x 6mm H x 308mm D
Interior light	N
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	SN
Refrigerant type	R134a
CFC free	Y
HFC free	N
HCFC free	Y
Ammonia free	Y
Optional plug types	Y
Voltage supply	230
Spark free	N

Auxiliary output	Y
Insulation thickness	40mm
Insulation type	Cyclopentane
Noise level	<45dB
Energy consumption	1.56 kw/24h
Door open alarm	Y
Air temperature alarm	Y
Antimicrobial handle	Y
Optional drawer system	N
Stacking kit compatible	N



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

