



## Under-counter Pharmacy Refrigerator Solid Door 107L PE307C

PE307C

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 2 year warranty as standard

### SPECIFICATION

Dimensions	495mm W x 830mm H x 620mm D
Weight	31.5kg
Product name	Pharmacy Fridge Solid Door 107L PE307C
GDPA code	PE307C
SKU code	UK 444441794, EU 444442337,
External dimensions	495mm W x 830mm H x 620mm D
Internal dimensions	403mm W x 700mm H x 400mm D
Packed dimensions	560mm W x 850mm H x 655mm D
Product weight	31.5kg
Packed weight	35kg
Gross capacity	107L
Net capacity	103L
External finish	Painted steel
Colour	White/green

Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door finish	Painted steel
Door style	Hinged reversible
Door open depth	1050mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	1 adjustable on front
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	2
Shelf material	Wire
Number of shelf positions	4
Shelf size	427mm W x 6mm H x 374mm D
Interior light	N
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	N
Refrigerant type	R134a
CFC free	Y
HFC free	N
HCFC free	Y
Ammonia free	Y
Optional plug types	Y
Voltage supply	230
Spark free	N

Auxiliary output	Y
Insulation thickness	35mm
Insulation type	Cyclopentane
Door open alarm	Y
Air temperature alarm	Y
Antimicrobial handle	Y
Optional drawer system available	N
Stacking kit compatible	N



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

