



## Countertop Pharmacy Refrigerator Glass Door 82L PG207C

PG207C

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 2 year warranty as standard

## SPECIFICATION

Dimensions	502mm W x 660mm H x 540mm D
Weight	30kg
Product name	Pharmacy Fridge Glass Door 82L PG207C
GDPA code	PG207C
SKU code	UK 444441793, EU 444442335,
External dimensions	502mm W x 660mm H x 540mm D
Internal dimensions	419mm W x 540mm H x 320mm D
Packed dimensions	555mm W x 677mm H x 580mm D
Product weight	30kg
Packed weight	33kg
Gross capacity	82L
Net capacity	67L
External finish	Painted steel
Colour	White/green

Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Glass
Door finish	Painted steel
Door style	Hinged right hand only, non reversible
Door open depth	976mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	N
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	2
Shelf material	Wire
Number of shelf positions	2
Shelf size	438mm W x 6mm H x 274mm D
Interior light	N
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	N
Refrigerant type	R134a
CFC free	Y
HFC free	N
HCFC free	Y
Ammonia free	Y
UK 3 pin plug	Y
Voltage supply	230
Spark free	N

Auxiliary output	N
Insulation thickness	30mm
Noise level	<43dB
Energy consumption	1.58 kw/24h
Door open alarm	Y
Air tempertaure alarm	Y
Antimicrobial handle	Y
Optional drawer system	N
Stacking kit compatible	N



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

