



Freestanding Large Pharmacy Refrigerator Glass Door 600L PGR600UK

PGR600UK

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display showing air and core temperatures
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ Temperature monitoring software records air and load temperatures constantly and also holds a full 7 day alarm history
- ✓ USB data download
- ✓ Power failure alarm

SPECIFICATION

Dimensions	690mm W x 1995mm H x 840mm D
Weight	140kg
Product name	Large Pharmacy Refrigerator Glass Door 600L PGR600UK
GDPA code	PGR600UK
SKU code	UK 444440365, EU 444442591
External dimensions	690mm W x 1995mm H x 840mm D
Internal dimensions	560mm W x 1230mm H x 680mm D
Packed dimensions	840mm W x 2195mm H x 990mm D
Product weight	140kg
Packed weight	155kg
Gross capacity	600L

External finish	Stainless steel
Colour	White/green
Internal finish	Stainless steel
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Glass
Door finish	Painted stainless steel
Door style	Hinged
Open door depth	1145mm
Replacement door seal	Y
Open door sensor	Y
Lock	Y
Adjustable feet	Fitted castors
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	5
Shelf material	Wire
Number of shelf positions	20
Shelf size	530mm W x 600mm D
Interior light	Y
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	ST
Refrigerant type	R134a
CFC free	Y
HFC free	Y
HCFC free	Y
Ammonia free	Y
UK 3 pin plug	Y

Voltage supply	230
Spark free	N
Auxiliary output	Y
Insulation thickness	60mm
Insulation type	Cyclopentane
Battery low alarm	Y (visual)
Air/load/defrost sensor fail alarm	Y (visual)
Compressor alarm	Y (visual)
Door open alarm	Y (visual and audio)
Air tempertaure alarm	Y (visual and audio)
Mains failure alarm	Y (visual and audio)
Antimicrobial handle	N
Optional drawer system available	N



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

