



## Freestanding Pharmacy Refrigerator Glass Door 273L PGRC273UK

PGRC273UK

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard

## SPECIFICATION

Dimensions	595mm W x 1565mm H x 625mm D
Weight	68.5kg
Product name	Pharmacy Fridge Glass Door 273L PGRC273UK
GDPA code	PSRC273UK
SKU code	444443778
External dimensions	595mm W x 1565mm H x 625mm D
Internal dimensions	480mm W x 1382mm H x 445mm D
Packed Dimensions	650mm W x 1585mm H x 727mm D
Product weight	68.5kg
Packed weight	73.5kg
Gross capacity	273L
Net capacity	269L
External finish	Painted aluminium/glass

Colour	White/green/glass
Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Painted aluminium/glass
Door finish	Painted aluminium/glass
Door style	Hinged, right handed only, non-reversible
Door open depth	1180mm
Replacement door seal	Y
Open door sensor	Y
Lock	Y
Adjustable feet	2 adjustable front feet, 2 rear rollers
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	5
Shelf material	Wire
Number of shelf positions	12
Shelf size	4 x 385mm W x 6mm H x 500mm D 1 x 340mm x 6mm H x 497mm D
Interior light	Y
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	SN
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
Ammonia free	Y
Optional plug types	Y

Voltage supply	220-240
Spark free	N
Auxiliary output	Y
Insulation thickness	60mm
Insulation type	Cyclopentane
Noise level	42
Energy consumption	1.08kW/24h
Door open alarm	Y
Air temperature alarm	Y
Antimicrobial handle	Y
Optional drawer system available	Y
Stacking kit compatible	N



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

