



Freestanding Pharmacy Refrigerator Glass Door 353L PGRC353UK

PGRC353UK

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard

SPECIFICATION

Dimensions	595mm W x 1860mm H x 640mm D
Weight	79kg
Product name	Pharmacy Fridge Glass Door 353L PGRC353UK
GDPA code	PSRC353UK
SKU code	444443786
External dimensions	595mm W x 1860mm H x 640mm D
Internal dimensions	480mm W x 1682mm H x 445mm D
Packed dimensions	650mm W X 1885mm H x 727mm D
Product weight	79kg
Packed weight	85kg
Gross capacity	353L
Net capacity	349L
External finish	Painted aluminium

Colour	White/green/glass
Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Painted aluminium/glass
Door finish	Painted aluminium/glass
Door style	Hinged right hand only, non reversible
Door open depth	1800mm
Replacement door seal	Y
Open door sensor	Y
Lock	Y
Adjustable feet	2 adjustable front feet, 2 rear rollers
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	6
Shelf material	Wire
Number of shelf positions	14
Shelf size	385mm W x 6mm H x 500mm D
Interior light	Y
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	SN
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
Ammonia free	Y
Optional plug types	Y
Voltage supply	220-240

Spark free	N
Auxiliary output	Y
Insulation thickness	60mm
Insulation type	Cyclopentane
Noise level	43dB
Energy consumption	1.39kW/24h
Door open alarm	Y
Air tempertaure alarm	Y
Antimicrobial handle	Y
Option drawer system available	Y
Stacking kit compatible	N



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

