



PSRC151UK Under-counter Pharmacy Refrigerator Solid Door 151L

PSRC151UK

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard

SPECIFICATION

Dimensions	595mm W x 845mm H x 600mm D
Weight	37kg
Product name	Pharmacy Fridge Solid Door 151L PSRC151UK
GDPA code	PSRC151UK
SKU code	444443767
External dimensions	595mm W x 845mm H x 600mm D
Internal Dimensions	500mm W x 705mm H x 405mm D
Packed dimensions	657mm W x 872mm H x 694mm D
Product weight	37kg
Packed Weight	41kg
Gross capacity	151L
Net capacity	141L
External finish	Painted steel
Colour	White/green
Internal finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door finish	Painted steel/solid
Door style	Hinged, reversible
Door open depth	1165mm
Replacement door seal	N
Open door sensor	Y
Lock	Y

Adjustable feet	1 adjustable front foot
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	3
Shelf material	Wire
Number of shelf positions	6
Shelf size	400mm W x 6mm H x 520mm D
Interior light	N
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	SN
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
Ammonia free	Y
Optional plug types	Y
Voltage supply	230
Spark free	N
Auxiliary output	Y
Insulation thickness	45mm
Insulation type	Cyclopentane
Noise level	41
Energy consumption	1.1kW/24h
Door open alarm	Y
Air temperature alarm	Y
Antimicrobial handle	Y
Optional draw system available	Y
Stacking kit compatible	Y



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

