



## Freestanding Pharmacy Refrigerator Solid Door 353L PSRC353UK

PSRC353UK

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard

## SPECIFICATION

|                     |   |
|---------------------|---|
| Dimensions          | 595mm W x 1860mm H x 640mm D              |
| Weight              | 69kg                                      |
| Product name        | Pharmacy Fridge Solid Door 353L PSRC353UK |
| GDPA code           | PSRC353UK                                 |
| SKU code            | 444443783                                 |
| External dimensions | 595mm W x 1860mm H x 640mm D              |
| Internal dimensions | 480mm W x 1682mm H x 445mm D              |
| Packed dimensions   | 650mm W x 1885mm H x 727mm D              |
| Product weight      | 69kg                                      |
| Packed weight       | 75kg                                      |
| Gross capacity      | 353L                                      |
| Net capacity        | 349L                                      |
| External finish     | Painted steel                             |

|                              |   |
|------------------------------|---|
| Colour                       | White/green                             |
| Internal finish              | HIPS*                                   |
| Type of plastic              | SRL800 (HIPS/ABS mix)                   |
| Door type                    | solid                                   |
| Door finish                  | Painted steel                           |
| Door style                   | Hinged, reversible                      |
| Door open depth              | 1185mm                                  |
| Replacement door seal        | N                                       |
| Lock                         | Y                                       |
| Adjustable feet              | 2 adjustable front feet, 2 rear rollers |
| Temperature range            | 2°C to 8°C                              |
| Temperature display on panel | Y                                       |
| Temperature alarm            | Y                                       |
| Number of shelves            | 6                                       |
| Shelf material               | Wire                                    |
| Number of shelf positions    | 14                                      |
| Shelf size                   | 385mm W x 6mm H x 500mm D               |
| Interior light               | Y                                       |
| Auto or manual defrost       | Auto                                    |
| Fan                          | Y                                       |
| Controller type              | Digital, external                       |
| Climate class                | SN                                      |
| Refrigerant type             | R600a                                   |
| CFC free                     | Y                                       |
| HFC free                     | Y                                       |
| HCFC free                    | Y                                       |
| Ammonia free                 | Y                                       |
| Optional plug types          | Y                                       |
| Voltage supply               | 220-240                                 |
| Spark free                   | N                                       |

|                         |              |
|-------------------------|--------------|
| Auxiliary output        | Y            |
| Insulation thickness    | 60mm         |
| Insulation type         | Cyclopentane |
| Noise Level             | 42dB         |
| Energy consumption      | 0.99kW/24h   |
| Door open alarm         | Y            |
| Air temperature alarm   | Y            |
| Antimicrobial handle    | Y            |
| Optional drawer system  | Y            |
| Stacking kit compatible | Y            |



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

**Wolflabs**

Tel: 01759 301142

Fax:01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

