



www.lec-medical.co.uk



## LC SRS236W TALL UPRIGHT LEC ESSENCHILL FRIDGE SOLID DOOR 236L WHITE

- ✓ **Large 236L Capacity**  
maximises available under-counter storage space
- ✓ **55cm Wide**  
For use where minimal space is needed
- ✓ **Energy Consumption 128 Kwh/year**
- ✓ **A+ Energy Rating**
- ✓ **Energy Efficient R600a Refrigerant**
- ✓ **Automatic Defrost**  
No need to manually defrost
- ✓ **3 Removable Glass Safety Shelves**
- ✓ **External Grab Handle**
- ✓ **2 Years Parts & Labour Warranty**



SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444444762

LC SRS236W

Tall Upright Lec essenChill Fridge Solid Door  
236L White



**FULL SPECIFICATION – LC SRS236W Tall Upright Lec essenChill Fridge  
Solid Door 236L Stainless Steel**

GDPA Code	LC SRS236W
SKU Number	444444762
External Product Dimensions (mm) (W x D x H)	550 x 550 x 1420
Product Weight (kg)	39
Packed Dimensions (mm) (W x D x H)	580 x 580 x 1462
Packed Weight (kg)	44
Gross Capacity (litres)	236
Net Capacity (litres)	235
Colour	White
Energy Rating	A+
Energy Consumption (kWh/annum)	128
Cost to Run (£/Year)	£20
Required Room Temperature	10 to 43 ° C
Climate Class	ST
Shelves (standard)	3 (safety Glass)
Door Storage	3
Salad Drawers	1
Refrigerant Type	R600a
Door Type	Solid
Removable Door Seals	Yes
Reversible Door	Yes
Handle Type (grab, integrated)	External Grab Handle
Self Closing Door	No
Defrost	Auto
Alarm Door Open	No
Noise Level	40 dB
Interior Light	Yes
Adjustable Thermostat	Yes
Levelling Feet	Yes
Warranty	2 Years Parts and Labour
Fitted UK 3 Pin Plug	Yes



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

