



## Labcold Flake Ice Makers



LITV-IQ50/IQ85C

**Labcold** Flake Ice Makers produce crushed ice ideal for filling packs and moulding to fit any item or object. The ice produced is both hygienic and clear so ideal for medical and laboratory applications.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

**Labcold** Ice Makers require both a water and electricity supply and a floor level drain.



LITV-IQ135

Part Number	LITV-IQ50	LITV-IQ85C	LITV-IQ135
Daily Ice Production	43kg per day	100kg per day	163kg per day
Ice Storage	14kg	20kg	60kg
Style	Freestanding	Freestanding	Freestanding
Exterior Dimensions (HxWxD)	750 x 400 x 520mm	800 x 470 x 600mm	1360 x 520 x 550mm
Interior Construction	Plastic Bins	Plastic Bins	Plastic Bins
Exterior Construction	Stainless Steel	Stainless Steel	Stainless Steel
Automatic Production	✓	✓	✓
Door Hinging	Bottom inward folding	Bottom inward folding	Bottom inward folding
Door Style	Solid	Solid	Solid
Feet	4	4	4



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

