

## RLHE1345 -40°C Sparkfree Chest Freezer



**Warning: Do not store open containers of volatile substances in this freezer.**

Robust chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

Part Number	RLHE1345
Capacity (Litres)	375
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1500 x 730mm
Interior Dimensions (HxWxD)	640 x 1295 x 445mm
Weight	68kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
Door Open Alarm	X
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Feet	4
Energy Consumption	0.8 kWh/24h
Eco Friendly Refrigerant	✓



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

