



RLVL03203 Basic Laboratory Freezer

Constructed from enamelled steel, this freezer from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** name and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability.

Part Number	RLVL03203
Capacity (Litres)	70
Style	Underbench
Exterior Dimensions (HxWxD)	835 x 500 x 490mm
Interior Dimensions (HxWxD)	370 x 305 x 365mm
Energy Consumption	0.54 kWh/24h
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -23°C
Sparkfree	✗
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Drawers	2 + 1 flap
Door Hinging	Right
Door Style	Solid
Feet	2 adjustable feet at front
Manual Defrost	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓



One clear fronted flap.



Two clear plastic drawers.

This 70 litre freezer has two clear plastic drawers and one flap located at the top. Two adjustable feet ensure it remains stable and manual defrost provides you with greater operator control over your temperature sensitive contents.





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

