



RLPL04043 Basic Laboratory Refrigerator

Part Number	RLPL04043
Capacity (Litres)	104
Style	Underbench
Exterior Dimensions (HxWxD)	841 x 495 x 500mm
Energy Consumption	0.323 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +10°C
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Shelves	3 + 1 bin
Sparkfree	✗
Door Hinging	Right
Door Style	Solid
Feet	2 adjustable feet at front
Internal Light	✓
Auto Defrost	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

Constructed from enamelled steel, this fridge from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** name and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability.



Two adjustable glass shelves.



One glass bin.

This 104 litre fridge has two adjustable shelves and one bin. It comes with an internal light and thermostat and has autodefrost for your convenience.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

