



## RLCF0720 Sparkfree Chest Freezer

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with a storage basket to assist with organisation and storage. Please ask at time of order for additional baskets. This chest freezer has manual defrost for greater operator control.

Part Number	RLCF0720
Capacity (Litres)	215
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 725 x 730mm (including Handle)
Interior Dimensions (HxWxD)	720 x 580 x 530mm
Weight	46kg
Energy Consumption	1.01kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	1
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
High Temperature Alarm	Visual
Manual Defrost	✓
Eco Friendly Refrigerant	✓





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

