



# LABCOLD™

medical & scientific refrigeration



## RLFF13248 Sparkfree Fridge Freezer

Part Number	RLFF13248
Capacity refrigerator	271 litres
Capacity freezer	113 litres
Style	Freestanding
Exterior Dimensions (HxWxD)	2000 x 595 x 650mm
Weight	67.5kg
Energy consumption	0.88kWh/24hr
Operating temperature	Fridge +0°C to +10°C Freezer -18°C to -23°C
Interior construction	Moulded Plastic
Exterior construction	White Enamel
Door style	Solid
Door hinging	Right
Door alarm	✓ Refrigerator only
Shelves	5
Drawers	3
Feet	4
Digital thermometer with temperature display, high and low alarms and min/max recording	Optional
High Temperature Alarm	✓ Visual only and only for the refrigerator compartment
Defrost	Fridge - auto Freezer - manual
Low GWP refrigerant	✓
1 year UK parts and labour warranty	✓

This Labcold Sparkfree fridge freezer is engineered to store flammable materials safely with all potential sources of ignition, such as the thermostat removed from inside the compartments of both the fridge and the freezer.

**Warning: Do not store open containers of volatile substances in this fridge freezer**

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board



Featuring a door open alarm on the refrigerator compartment as standard, this fridge freezer can also be factory fitted with a Labcold digital thermometer which has a high and low alarm and records minimum and maximum temperatures. If you require door locks, these can also be fitted at time of order, simply use RLFF13248LK as the product code.

As you would expect, your sparkfree fridge freezer is backed by a comprehensive UK parts and labour warranty and as Labcold are an accredited UKAS temperature calibration laboratory #8898 we can take care of all your calibration needs as part of a general maintenance contract. Please ask for a quote when you place your order.





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

