



LABCOLD™

medical & scientific refrigeration



RLPR1514 Sparkfree Fridge

Part Number	RLPR1514
Capacity (Litres)	400
Style	Free standing
Exterior Dimensions (HxWxD)	1855 x 600 x 585mm
Interior Dimensions (HxWxD)	1620 x 506 x 485mm
Weight	69kg
Energy Consumption	1.34kWh/24hr
Operating Temperature	0°C to +10°C
Factory Set Point	+4°C
Interior Construction	Moulded ABS
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right (Can be reversed at time of manufacture)
Door Lock	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	5 at 500 x 410mm 1 at 500 x 230mm
Feet	2 (adjustable)
Rollers	2 at rear
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Min/Max Recording	✓
Access Port	Optional Extra
Automatic Defrost	✓
Fan Assisted Cooling	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm. The controller also records the minimum and maximum temperature since the last reset.





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

