

## ROTARY-MIXER

**MODEL - L26** 50/60Hz CAT.NO.: 26000 (230V) / 26010 (115V)

- Provides serial rotation of test-tubes
- With diameter-range: 8 to 27 mm
- All aluminium and stainless-steel construction
- Speed: stepless controllable: 5 to 50 r.p.m.
- Quick coupling / exchange of rotating drum / bar
- Provided with stainless-steel drip-tray.
- Delivered with bar for 20 tubes (14-18 mm)
- Dimensions: 540x160x185 mm Weight: 7,5 Kg.



### ACCESSORIES FOR ROTARY-MIXERS

Bar/Drum outfit	tube Ø range mm	cat.no. for L26 and L27	max.tube length mm		max no. tubes
			L-26	L-27	
	14-18	26020	190	290	20
	18-22	26022	185	290	18
	22-26	26027	180	290	18
	8-12	26040	85	135	104
	11-15	26480	120	180	48
	If 2 rows used		140	250	24
	8-12	26560	120	180	56
	If 2 rows used		140	250	28

### SPECIAL / LARGE ACCESSORIES

	CATALOG NO: 27304 BOTTLE ROLLING ACCESSORY For rolling two or four bottles of max length 300 mm Diameter range 25 to 120 mm
	CATALOG NO: 27312 HYBRIDISATION ACCESSORY For rolling maximum 12 hybridisation-bottles Diam: 35 mm Length 300 mm
	CATALOG NO: 27400 BOTTLE ACCESSORY FOR TOP-OVER-BOTTOM-TURNING of bottles. Diam. range: 25-125 mm Length range: 60-170 mm.
	CATALOG NO: 27580 EXTRACTION ACCESSORY For measuring flasks. Max. 18 flasks 50ml.

## LARGE - ROTARY - MIXER

**MODEL-L27** 50/60Hz CAT.NO.: 27000 (230V) / 27010 (115V)

- Same as MODEL L26 but 50 mm higher
- Provides rotation action for up to 310 mm diameter accessories





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

