

TEST - TUBE - ROTATORS



VARIABLE-SPEED MODEL - L29

50/60Hz CAT.NO.: 29000 (230V) / 29010 (115V)

- Same as L-28 BUT VARIABLE SPEED 5 TO 50 R.P.M.
- Dimensions: 210(l) x 200(d) x 270 (h)
- Weight: 3,6 Kg. (without disk)



ROTATING DISKS

catalog number	number clamps	tube Ø range
28012	20	8-12
28020	20	11-15
28030	20	14-18
28040	10	32-38
28050	15	22-26
28060	40	4-6
28070	18	18-22
28080	12	26-32
28090	adapter for dual disk support	



MODEL - L28

50/60Hz CAT.NO.: 28000 (230V) / 28010 (115V)

- For mixing, diffusion, dialysing, extraction etc.
- SPEED 10 r.p.m. With strong, geared motor.
- Angle of rotation can be set between 90° and 27°
- Tubes may be exchanged while rotator is running
- Comes with line-cord, but without disks.
- Dimensions: 210(l) x 200(d) x 220 (h), Weight: 2,6 kg.



ROCK 'N ROLLER MODEL - L201

50/60Hz CAT.NO.: 20100 (230V) / 20110 (115V)

- Keeps samples in gentral motion by rocking while rolling.
- 14 stainless steel rotating bars, 208 x 10 mm diam.
- with a heart to heart distance of 19,5 mm
- Speed of rotation is approx. 10 r.p.m.
- Rocking frequency is approx. 5 cycles / min. at 35° movement.
- Front panel mounted, illuminated on-off switch.
- Dimensions: 330(l) x 300(d) x 190 (h).
- Weight: 3,0 kg.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

