


MODEL LD-76

ROTARY MIXER




- Provides serial rotation of tubes and mini-tubes for extraction, diffusion, etc.
- Mixing action may be set by speed control 5-50 rpm.
- Tube-diameter-range 4-32 mm. Chosen from listing below.
- Quick/easy coupling/exchange of rotation drums/bars.
- Tubes may be exchanged while the rotator is running.
- Dimensions: LxWxH 450x170x190 mm.
- Weight: 3,3 Kg.

50/60 HZ CAT.NO. 76000 (230V) / 76010 (115V)

ELECTRONICALLY CONTROLLED STEPPER-MOTOR PROVIDES:

- Digital read-out for - speed 5-50 rpm.
- time 0-99 hrs/1-59 min.
- Electronic timer for "auto-stop".
- "Memory-function" for repeated operation.
- "Re-start-function" for continued operation after line-voltage failure including auto-stop.
- "Settings saving" available.

- ON/OFF MAINS
- START/STOP MOTOR
- "MODE" SWITCH
- "UP/DOWN"
- "TIMER/CYCLING"



ACCESSOIRES FOR ROTARY-MIXERS

Bar/Drum outfit	tube \varnothing range mm	cat.no. for L26/LD-76	max.tube length mm	max no. tubes
	16-20	26020	190	20
	19-22	26022	185	18
	22-27	26027	180	18
	4-6	26036	85	104
	6-9	26039	85	104
	9-12	26040	85	104
	12-16	26480	120	48
	if 2 rows used		140	24
	9-12	26560	120	56
	if 2 rows used		140	28



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

