



50/60 HZ CAT.NO. 79000 (230V) / 79010 (115V)

## ELECTRONICALLY CONTROLLED STEPPER-MOTOR PROVIDES:

- Digital read-out for - speed 5-50 rpm.  
- time 0-99 hrs/1-59 min.
- Electronic timer for "auto-stop".
- "Memory-function" for repeated operation.
- "Re-start-function" for continued operation after line-voltage failure including auto-stop.
- "Settings saving" available.




## MODEL L-28 TEST-TUBE ROTATOR

- Provided with strong 10 rpm. synchronous motor
- Accepts all disks mentioned in the listing besides
- Dimensions (LxWxH) 210x200x220 mm.
- Weight: 3,5 Kg.

50/60 HZ CAT.NO. 28000 (230V) / 28010 (115V)

## MODEL LD-79 TEST-TUBE-ROTATOR

- Provides mixing for diffusion, dialysing, extraction etc.
- Tube diameter range 4-40 mm.
- Accepts all disks mentioned in the listing below.
- Mixing action can be set by:
  - speed control 5-50 rpm.
  - angle of rotation 90° to 27°.
- Quick/easy exchange of disks by unscrewing central nut.
- Dual or triple disk support may be carried out.
- Tubes may be exchanged while the rotator is running.
- Dimensions: 200 (L) x 200 (W) x 250 (H) mm.
- Weight: 3.5 Kg.



ROTATING DISKS				FOR TEST TUBES		
FOR MINI-REACTION VIALS/TUBES				CATALOG NUMBER	NUMBER CLAMPS	TUBE Ø RANGE
CATALOG NUMBER	NUMBER CLAMPS	TUBE Ø RANGE	VOLUME ML			
28006	40	4-6	0.4	28020	20	12-16
28009	30	6-9	0.5	28030	20	16-19
28012	20	9-12	1.5	28032	10	28-32
28090	adapter for dual disk support			28040	10	32-40
				28050	15	22-25
				28070	17	19-22
				28080	12	25-28



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

