

Thermo Mix

HC110-Pro

Dry bath with cooling

Features:

- Precise temperature control for heating and cooling
- Flexible adapter selection
- Programmable
- Compatible with various tubes



Excellent temperature accuracy



Programmable



Adapter self recognition



Cat.No.	Blocks	Speed	Capacity	Diameter x Deep
18900238		1500rpm	24 x 0.5ml	8.1x26mm
18900237		1500rpm	24x 1.5ml	10.9x30.8mm
18900239		1500rpm	24x 2ml	10.9x30.8mm
18900240		1500rpm	8x 5ml	13.6x42mm
18900241		800rpm	8x 15ml	16.5x80mm
18900242		600rpm	4x 50ml	29.6x80mm
18900244		1500rpm	well plate	-
18900243		1500rpm	0.2ml x 96 PCR plate	-

Specifications

	HC110-Pro
Functions	Heating& cooling
Temperature range	Heating: room temp. - 110°C Cooling: room temp. - below room temp. 25°C
Temperature control accuracy[@ 20-45°C]	± 0.5°C
Temperature uniformity[@ 20-45°C]	Max ± 0.5 °C
Max heating rate	5.5°C/min
Max cooling rate	2.5°C/min (100°C-room temp.) 0.5°C/min (below room temp.)
Screen	TFT
Program	9
Overheating protection	150°C
Adapter block material	aluminum
Voltage, Frequency	100-240V,50/60Hz
Power	200W
Dimension[WxDxH] (without the heating block)	200x235x120mm
Weight	7.3kg



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

