



Countertop Laboratory Spark Free Refrigerator 39L

GDPA code: **LSFSF39UK**

SKU code: **444410574**

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



The Experts in
Medical Refrigeration

Specification: Countertop Laboratory Spark Free Refrigerator 39L

Gross Capacity (litres)	39
Net Capacity (litres)	35
Unit External Width (mm)	480
Unit External Depth (mm)	505
Unit External Height (mm)	565
Product Weight (kg)	18
Packed Width (mm)	490
Packed Depth (mm)	530
Packed Height (mm)	620
Packed Weight (kg)	21
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	360 x 340 x 395
Temperature Range	-18°C to -25°C (tolerance + or - 2)
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	905
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	4 - All adjustable
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear LH Corner
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	1
Number of Shelf Positions	1
Shelf Size (mm) (WxDxH)	365 X 190
Shelf Material	Glass
Number of Baskets	0
Basket Material	N/A
Basket Size (mm) (WxDxH)	363 x 150
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Auto or Manual Defrost	None
Controller Type	Digital, External
Number of Probes	1
Probe Type	Air
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes

Specification: Countertop Laboratory Spark Free Refrigerator 39L

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	—
Energy Consumption (kwh)	—
Voltage Supply	230v
Aux Output	No
Insulation Thickness (mm)	50
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	Yes



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

