



## Under-counter Laboratory Spark Free Freezer 98L

GDPA code: **LSFSF98UK**

SKU code: **444410535**

- ✓ -18°C to -25°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms




The Experts in  
Medical Refrigeration

## Specification: Under-counter Laboratory Spark Free Freezer 98L

<b>Gross Capacity (litres)</b>	98
<b>Net Capacity (litres)</b>	84
<b>Unit External Width (mm)</b>	595
<b>Unit External Depth (mm)</b>	642
<b>Unit External Height (mm)</b>	850
<b>Product Weight (kg)</b>	37.25
<b>Packed Width (mm)</b>	640
<b>Packed Depth (mm)</b>	660
<b>Packed Height (mm)</b>	910
<b>Packed Weight (kg)</b>	41.75
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	420 x 430 x 640
<b>Temperature Range</b>	-18°C to -25°C (set point tolerance + or - 2)
<b>Colour</b>	White/Green
<b>External Finish</b>	Painted Steel
<b>Internal Finish</b>	HIPS*
<b>Type of plastic</b>	SRL800 (HIPS/ABS mix)
<b>Door Finish</b>	Painted Steel / Solid
<b>Door Type</b>	Solid
<b>Door Style</b>	Hinged, Reversible
<b>Door Open Depth (mm)</b>	1130
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	Yes
<b>Lock</b>	Yes, 2 Keys
<b>Adjustable Feet (number &amp; position)</b>	4 - All adjustable
<b>Porthole</b>	Yes
<b>Porthole Size</b>	16 mm
<b>Porthole Location</b>	Rear LH Corner
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	Yes
<b>Remote Alarm Contacts</b>	Yes
<b>Audible Alarms</b>	Yes
<b>Visual Alarms</b>	Yes
<b>Number of Shelves</b>	3
<b>Number of Shelf Positions</b>	3
<b>Shelf Size (mm) (WxDxH)</b>	495 x 325
<b>Shelf Material</b>	Glass / Plastic
<b>Number of Baskets</b>	0
<b>Basket Material</b>	N/A
<b>Basket Size (mm) (WxDxH)</b>	485 x 180 x 171
<b>Interior Light</b>	No
<b>Light Type</b>	N/A
<b>Light Activation</b>	N/A
<b>Fan</b>	No
<b>Spark Free</b>	Yes
<b>Auto or Manual Defrost</b>	None
<b>Controller Type</b>	Digital, External
<b>Number of Probes</b>	1
<b>Probe Type</b>	Air
<b>Climate Class</b>	SN
<b>Refrigerant Type</b>	R600a

**Specification:** Under-counter Laboratory Spark Free Freezer 98L



<b>CFC Free</b>	Yes
<b>HFC Free</b>	Yes
<b>HCFC Free</b>	Yes
<b>Ammonia Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	230v
<b>Insulation Thickness (mm)</b>	85
<b>Insulation Type</b>	High Density Polyurothane Foam
<b>Blowing Agent</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes
<b>Optional Drawer System Available</b>	No
<b>Wall Brackets</b>	No
<b>Digital Lock</b>	No



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

