



Under-Counter Laboratory Control Plus Fridge-Freezer Combi 119L ATEX

GDPA code: LSFC119UK-ATEX Lab ATEX UC F/F

SKU code: 444410685

- ✓ Operating Temperature Range of +2°C to +10°C for Fridge and -18°C to -25°C for Freezer
- ✓ Spark Free
- ✓ Advanced Controller
- ✓ Unique Anti Microbial Handle
- ✓ Fitted Lock as Standard



The Experts in
Medical Refrigeration

Specification: Under-Counter Laboratory Control Plus Fridge-Freezer Combi 119L ATEX

Gross Capacity (litres)	119 (Fridge: 106, Freezer: 13)
Net Capacity (litres)	117L (fridge 105L, freezer 12L)
Unit External Width (mm)	540
Unit External Depth (mm)	600
Unit External Height (mm)	860
Product Weight (kg)	34
Packed Width (mm)	602
Packed Depth (mm)	639
Packed Height (mm)	877
Packed Weight (kg)	36
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	Fridge: 450 x 450 x 525, Freezer: 400 x 295 x 120 (step: 155d x 225h)
Temperature Range	Fridge: +2°C to +10°C, Freezer: -18°C to -25°C
Colour	White/Green
External Finish	White
Internal Finish	White
Type of plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1090
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet (number & position)	2 Adjustable, 2 Rear Rollers
Porthole	-
Porthole Size	-
Porthole Location	-
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	2
Number of Shelf Positions	4
Shelf Size (mm) (WxDxH)	308 x 448 x 4.5, 1: 200 x 448 x 4.5
Shelf Material	Wire
Number of Baskets	0
Basket Material	1 Solid
Basket Size (mm) (WxDxH)	2 full + 2 half
Interior Light	N/A
Light Type	No
Light Activation	-
Fan	-
Spark Free	No
Auto or Manual Defrost	Yes
Controller Type	Fridge: Automatic, Freezer: Manual
Number of Probes	Digital, External
Probe Type	1
Climate Class	Air
Refrigerant Type	N
CFC Free	R600a
HFC Free	Yes

Specification: Under-Counter Laboratory Control Plus Fridge-Freezer Combi 119L ATEX

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	Yes
Energy Consumption (kwh)	-
Voltage Supply	-
Aux Output	230
Insulation Thickness (mm)	Yes
Insulation Type	45
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

