



British healthcare and laboratory
refrigeration specialists



Laboratory Range

Under Counter Fridge and Freezer



FEATURES

Spark free interior

Lock

12 month warranty

KEY BENEFITS

Our range of laboratory refrigerators and freezers offers the safe and reliable storage of volatile liquids, reagents, chemicals, poisons and samples with a spark-free interior that reduces the risk of internal explosion.

Security is essential in laboratories so our Lec Medical fridge and freezer come complete with door lock and two keys, ensuring that only authorised individuals are able to access them.

Our under counter fridge and freezer come complete with a 12 month warranty for added peace of mind.




British healthcare and laboratory
refrigeration specialists

Laboratory Range

Under Counter Fridge and Freezer



SPECIFICATION TABLE



Product Description		Laboratory Fridge
Product Dimensions (external)	mm	545 x 554 x 845mm (W x D x H)
Capacity	litres	130
Colour/Finish		White
Reversible Door	Y/N	Yes
Door Lock	Y/N	Yes
Drawer Number		1 Plastic
Shelf Number		3 Glass
Handle Type		Recessed
Defrost Capability		Auto Defrost
Annual Energy Consumption	kwh	117
Thermostat		Adjustable
Spark Free	Y/N	Yes
HFC Free	Y/N	Yes
CFC Free	Y/N	Yes
Product Description		Laboratory Freezer
Product Dimensions (external)	mm	545 x 554 x 845mm (W x D x H)
Capacity	litres	85
Colour/Finish		White
Reversible Door	Y/N	Yes
Door Lock	Y/N	Yes
Drawer Number		3 Plastic
Shelf Number		1 Plastic
Handle Type		Recessed
Annual Energy Consumption	kwh	157
Thermostat		Adjustable
Spark Free	Y/N	Yes
HFC Free	Y/N	Yes
CFC Free	Y/N	Yes



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

