



Countertop Laboratory Essential Refrigerator 82L

GDPA code: **LSFSR82UK**

SKU code: **444410643**

- ✓ 2°C to 10°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



The Experts in
Medical Refrigeration

Specification: Countertop Laboratory Essential Refrigerator 82L

Gross Capacity (litres)	82
Net Capacity (litres)	67
Unit External Width (mm)	502
Unit External Depth (mm)	660
Unit External Height (mm)	540
Product Weight (kg)	27.5kg
Packed Width (mm)	555
Packed Depth (mm)	677
Packed Height (mm)	580
Packed Weight (kg)	30.5kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	419mm W x 540mm H x 315mm D (215mm H x 165mm D)
Temperature Range	2°C to 10°C
Colour	White/green
External Finish	White painted steel
Internal Finish	White plastic
Type of plastic	HIPS/ABS mix
Door Finish	Painted steel
Door Type	Solid
Door Style	Hinged, reversible
Door Open Depth (mm)	976mm
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	No
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	2
Number of Shelf Positions	2
Shelf Size (mm) (WxDxH)	438mm W x 6mm H x 274mm D
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Basket Size (mm) (WxDxH)	—
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Auto or Manual Defrost	Auto
Controller Type	Digital, External
Number of Probes	1
Probe Type	Air
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	No

Specification: Countertop Laboratory Essential Refrigerator 82L

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	<43dB
Energy Consumption (kwh)	0.77 kW/24h
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	30mm
Insulation Type	Cyclopentane
Blowing Agent	Yes
Anti Microbial Handle	No
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

