



Advanced Temperature Probe Data Logger

ATMDL-LCDGP

- ✓ -40 to +125°C (-40 to +257°F) measurement range
- ✓ Stores over 32,000 readings
- ✓ Logging rates between 1 second and 12 hours
- ✓ High contrast LCD, with four digit temperature display function
- ✓ Immediate, delayed and push-to-start logging
- ✓ User-programmable alarm thresholds
- ✓ Status indication via red and green LEDs
- ✓ Supplied with a Type 2 thermistor probe

SPECIFICATION

Dimensions	135mm W x 24mm H x 21mm D
Weight	N/A
Measurement range	-40°C to 125°C
Temperature display on panel	Y
Accuracy	+/-0.1°C



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

