

## Digital Mantles (DM Series)

**Brand:** Medline Scientific

The DM-Series feature an integral digital temperature controller linked to a sensor probe, which is inserted directly into the flask. The temperature is shown clearly on an LED display and the desired temperature is controlled by a simple push button potentiometer. A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



## Specifications:

Product Code(s)	DM-601	DM-602	DM-603	DM-604	DM-605	DM-606	DM-607	DM-608	DM-609
Case Construction	Steel, with chemical resistant coating								
Thermal Insulation	Moulded silica powder								
Heating Element	Nichrome wire								
Mantle	Woven glass fibre mat								
Element Temperature	Up to 450°C (max.)								
Sensor Probe	CA (K) Type, Ø 4.8mm * 200mm (1.5m lead)								
Temperature Range (options)	Low (up to 100°C)/High (up to 200°C)								
Power Supply	220V AC								
Volume (ml)	100	250	500	1000	2000	3000	5000	10,000	20,000
Watts	80 W	150 W	260 W	350 W	450 W	550 W	700 W	1200 W	1800 W

ISO 9001 CERTIFIED ORGANISATION



United Registrar of Systems Certificate No. 1938



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

